

## Me-Hot Thermal Imaging Technical Training



Trainings are held at the Med-Hot Training Facility in Lakeland, FL., conveniently located in the center of the state between Orlando and Tampa.

This is the basic technical user training offered to purchasers of the Med-Hot Systems. We teach our users how to use the TotalVision Software and how to take all images to internationally accepted protocols. We can effectively train two hands-on technicians but others from your office may come along to watch and learn.

Trainings are scheduled for 4 days....ending at noon on the fourth day to allow attendees time for travel later that day.

Healthy lunches and snacks will be provided on site. Please let us know of any special food sensitivities or preferences.

### Preparation:

We will provide a:

- Complete User Manual

- Thirty minute Introduction Video to be reviewed before the course

- List of supplies for the clinic

Client provides a:

- List of each person's name what each person's role will be in the imaging process. Identify the primary thermographer.

- Travel arrival and departure information and attendees cell phone

Subjects covered:

- How to set up the clinic for privacy imaging

- Understanding the physiology of the skin and it's functional association with the various systems of the body

- How to represent thermography to the patient and in advertising

- Patient and environmental protocols

- Review TotalVision Software and hardware functionality

- Imaging protocols for each region of interest

Hands-on Training – Multiple models will be available for learning and practice

- Demonstration

- Using the software and features

- Practice camera movement and the patient

- Practice focusing and zoom

- Step-by-step review and practice of each region

- Each trainee will be responsible for at least one completed full body study

# **Me-Hot Thermal Imaging Technical Training**

## **Review and Testing**

We provide a test review on the 3<sup>rd</sup> day

Review studies and self-critique

Practical test consisting of a complete full body with no assistance from the trainer

### **Written Test**

Each trainee will be assigned a number of studies to submit after returning home. It is required that trainees have arranged for patients, friends or family to be available for practice immediately after returning home.

A Certificate of Training will be provided