

# Med-Hot Automated Stand

## Essential for Privacy Health Screening

### Professional, Easy-To-Use, Affordable

**med-hot**

Health Discovery with Infrared

Specially engineered for ease of use and privacy screening. Available as a portable mobile stand or fixed in-office wall mount. Conveniently control vertical camera positioning from the desktop through an easy to use on-screen graphical user interface.

Key Features	Benefits
<ul style="list-style-type: none"><li>• Computer On-screen Control</li></ul>	Efficient remote control eliminates manual adjustment. Easily fine-tune imaging camera position.
<ul style="list-style-type: none"><li>• Pan-Tilt Positioning</li></ul>	Reduce session imaging time.
<ul style="list-style-type: none"><li>• Prewired</li></ul>	Quick set-up and tear down.
<ul style="list-style-type: none"><li>• Camera Quick Release</li></ul>	Quick set-up and tear down. Quickly install and de-install the camera & pan-tilt motor on the camera platform.
<ul style="list-style-type: none"><li>• Portable (Mobile Stand)</li></ul>	Compact and easy to transport. Three wheel collapsible stand easily folds up for travel.



Portable Mobile Stand



In-office Wall Mount


# Med-Hot Automated Stand

## Essential for Privacy Health Screening

### Professional, Easy-To-Use, Affordable

# med-hot

Health Discovery with Infrared



Desktop On-screen Graphical User Interface

## Specifications

	<b>Automated Mobile Stand</b>	<b>Wall Mount</b>
• Actuator Stroke	30 in (76.2 cm)	40 in (101.6 cm)
• Maximum Camera Height	48.5 in (123.2 cm)	58.5 in (148.59 cm) Note: Can vary dependent upon installation height.
• Minimum Camera Height	18 in (45.7 cm)	18 in (45.7 cm) Note: Can vary dependent upon installation height.
• Overall Height	55 in (139.7 cm) Unfolded 51 in (129.5 cm) Folded	51 in (129.5 cm) Note: Can vary dependent upon installation height.
• Leveling Leg (1)	Yes	N/A
• Weight	27 lbs (12.2 kg)	
• Input Power	100- 240 V, XX A, 50-60 Hz	100- 240 V, XX A, 50-60 Hz

### Contact Us

For more information, visit [Med-Hot.com](http://Med-Hot.com) or call 1 800.646.1599 / 1 800.709.9565 to speak to a representative.

Copyright 2017 Med-Hot Thermal Imaging, In. All rights reserved. Med-Hot is a registered trademark of Med-Hot Thermal Imaging, Inc. The names and logos of any companies product mentioned herein may be the trademarks of their respective owners in the United States, Canada, and/or other countries. The information herein is subject to change without notice.