

Med-Hot Portable Manual Stand

Basic Portable Stand, Manual Positioning,
Economical, Entry-Level

med-hot

Health Discovery with Infrared

Basic portable stand with manual positioning is an excellent economical option for entry-level thermographers. Includes fold up legs and manual pan-tilt platform.

Key Features	Benefits
<ul style="list-style-type: none">• Braked wheels	Stable base. Three wheels lock to prevent movement.
<ul style="list-style-type: none">• Fold-up Legs	Easy to set-up and break down for transport
<ul style="list-style-type: none">• Manual Pan-Tilt with Adjustable Handle	Ease of use. Smooth balanced movement with locks.
<ul style="list-style-type: none">• Spring Loaded Sliding Arm	Quick height adjustment.



Portable Manual Stand



Manual Pan-Tilt

Med-Hot Automated Stand

Basic Portable Stand, Manual Positioning,
Economical, Entry-Level

med-hot

Health Discovery with Infrared

Specifications

• Footprint diameter	30 in (76.2 cm)
• Maximum camera height	52.5 in (133.4 cm)
• Minimum camera height	14.5 in (36.8 cm)
• Overall height	52 in (132.1 cm)
• Collapsed height	48 in (121.9 cm)
• Weight	14.25 lbs (6.46 kg)

Contact Us

For more information, visit Med-Hot.com or call 1 800.646.1599 / 1 800.709.9565 to speak to a representative.

Copyright 2017 Med-Hot Thermal Imaging, In. All rights reserved. Med-Hot is a registered trademark of Med-Hot Thermal Imaging, Inc. The names and logos of any companies product mentioned herein may be the trademarks of their respective owners in the United States, Canada, and/or other countries. The information herein is subject to change without notice.